

Product datasheet for **TP726565**

TCFL5 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transcription factor-like 5 protein/TCFL5
Species:	Human
Expression cDNA Clone or AA Sequence:	A139-L425
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Transcription factor-like 5 protein/TCFL5is produced by E. coli. The target gene encodingA139-L425 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10732
UniProt ID:	Q9UL49
Synonyms:	Transcription factor-like 5 protein; Cha transcription factor; HPV-16 E2-binding protein 1 (E2BP-1); TCFL5; CHA; E2BP1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)