

Product datasheet for **TP726568**

ZNF763 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Zinc finger protein 763/ZNF763
Species:	Human
Expression cDNA Clone or AA Sequence:	S268-L394
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Zinc finger protein 763/ZNF763is produced by E. coli. The target gene encodingS268-L394 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	284390
UniProt ID:	Q0D2J5



[View online »](#)