

Product datasheet for **TP726569**

RPS7 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Small ribosomal subunit protein eS7/RPS7
Species:	Human
Expression cDNA Clone or AA Sequence:	F2-L194
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Small ribosomal subunit protein eS7/RPS7 is produced by E. coli. The target gene encodingF2-L194 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	6201
UniProt ID:	P62081



[View online »](#)