

Product datasheet for **TP726570**

MTHFS Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human 5-formyltetrahydrofolate cyclo-ligase/MTHFS
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-A203
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human 5-formyltetrahydrofolate cyclo-ligase/MTHFS is produced by E. coli. The target gene encoding A2-A203 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10588
UniProt ID:	P49914



[View online »](#)