

Product datasheet for **TP726571**

PGRMC1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PGRMC1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	K44-D195
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human PGRMC1 Protein is produced by E. coli. The target gene encoding K44-D195 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10857
UniProt ID:	<u>O00264</u>



[View online »](#)