

Product datasheet for **TP726571**

PGRMC1 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human PGRMC1 Protein |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | K44-D195 |
| Tag: | C-FLAG |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human PGRMC1 Protein is produced by E. coli. The target gene encoding K44-D195 is expressed with a FLAG tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 10857 |
| UniProt ID: | <u>O00264</u> |



[View online »](#)