

Product datasheet for **TP726572**

TMEM40 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Transmembrane protein 40/TMEM40
Species:	Human
Expression cDNA Clone or AA Sequence:	E2-E159
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Transmembrane protein 40/TMEM40is produced by E. coli. The target gene encodingE2-E159 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	55287
UniProt ID:	<u>Q8WWA1</u>



[View online »](#)