

## Product datasheet for **TP726573**

### p21 ARC (ARPC3) Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human Actin-related protein 2/3 complex subunit 3/ARPC3  |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | P2-Q178  |
| Tag:                                  | C-FLAG   |
| Buffer:                               | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | Recombinant Human Actin-related protein 2/3 complex subunit 3/ARPC3is produced by E. coli. The target gene encodingP2-Q178 is expressed with a FLAG tag at the C-terminus. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 10094  |
| UniProt ID:                           | <u><a href="#">O15145</a></u>  |



[View online »](#)