

## Product datasheet for **TP726576**

### **PITPNB Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant Human PITPNB Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	l111-V271
Tag:	C-FLAG
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human PITPNB Protein is produced by E. coli. The target gene encoding l111-V271 is expressed with a FLAG tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	23760
UniProt ID:	<a href="#">P48739</a>



[View online »](#)