

Product datasheet for **TP726577**

FOXS1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Forkhead box protein S1/FOXS1
Species:	Human
Expression cDNA Clone or AA Sequence:	H57-G303
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Forkhead box protein S1/FOXS1 is produced by E. coli. The target gene encoding H57-G303 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2307
UniProt ID:	O43638
Synonyms:	Forkhead box protein S1; Forkhead-like 18 protein; Forkhead-related transcription factor 10 (FREAC-10); FOXS1; FKHL18; FREAC10
Protein Families:	Transcription Factors



[View online »](#)