

## Product datasheet for **TP726578**

### POU5F2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SPRM-1/POU5F2
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-Q144
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human SPRM-1/POU5F2is produced by E. coli. The target gene encodingA2-Q144 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	134187
UniProt ID:	<a href="#">Q8N7G0</a>
Synonyms:	POU domain, class 5, transcription factor 2; Sperm 1 POU domain transcription factor (SPRM-1); POU5F2; SPRM1
Protein Families:	Transcription Factors



[View online »](#)