

Product datasheet for **TP726585**

CD3E Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CD3E (C-mFc)
Species:	Human
Expression cDNA Clone or AA Sequence:	Gln22-Thr48
Tag:	C-mFc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Note:	Recombinant Human T-Cell Surface Glycoprotein CD3 epsilon Chain is produced by our Mammalian expression system and the target gene encoding Gln22-Thr48 is expressed with a mFc tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	916
UniProt ID:	P07766
Summary:	T-Cell Surface Glycoprotein CD3 ϵ Chain (CD3 ϵ) is a single-pass type I membrane protein. CD3 ϵ contains 1 Ig-like (immunoglobulin-like) domain and 1 ITAM domain. CD3 ϵ is a polypeptide encoded by the CD3E gene on chromosome 11 in humans. The T cell receptor-CD3 complex (TCR/CD3 complex) is involved in T-cell development and several intracellular signal-transduction pathways. This complex is critical for T-cell development and function, and represents one of the most complex transmembrane receptors. The T cell receptor-CD3 complex is unique in having ten cytoplasmic immunoreceptor tyrosine-based activation motifs (ITAMs). TCR/CD3 complex plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways.



[View online »](#)