

Product datasheet for TP726652

TSLP Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant Human TSLP (R127A, R130A, C-10His)

Species: Human

Expression cDNA Clone or AA Sequence: Tyr29-Gln159(R127A, R130A)

Tag: C-10His

Buffer: Lyophilized from a 0.2 um filtered solution of PBS, pH7.4.

Note: Recombinant Human Thymic stromal lymphopoietin is produced by our Mammalian expression system and the target gene encoding Tyr29-Gln159(R127A, R130A) is expressed with a 10His tag at the C-terminus.

Stability: 12 months from date of despatch

Locus ID: 85480

UniProt ID: [Q969D9](#)

Summary: Thymic stromal lymphopoietin (TSLP) is a novel member of the hemopoietic cytokine family that promotes the development of B cells and shares overlapping activity with IL-7. The human TSLP protein comprises a 28 amino acids (aa) signal sequence and 131 aa mature region. Human TSLP has two isoforms lftTSLP and sfTSLP produced by alternative splicing . lftTSLP is expressed in a number of tissues including heart, liver and prostate, and sfTSLP (63aa) is predominantly expressed in keratinocytes of oral mucosa, skin and in salivary glands. In aa sequence level, Human TSLP displays about 43% identity with mouse TSLP. TSLP is a cytokine that functions mainly on myeloid cells; it induces the release of T cell-attracting chemokines from monocytes and enhances the maturation of CD11c(+) dendritic cells. TSLP has proliferative effects on the myeloid cell line and may initiate asthma or atopic dermatitis responses by directly activating mast cells . TSLP signals cells via the interleukin-7 receptor- α chain (IL-7R α), shared with IL-7, together with the TSLP receptor (TSLPR) subunit. Recent studies indicate that TSLP and its receptor are novel therapeutic targets for rheumatoid arthritis, for increased intraarticular TSLP concentrations in patients has caused chemotaxis and activation of arthritogenic T cells.



[View online »](#)