

Product datasheet for **TP726669**

Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human IL-13 (C-mFc)
Species:	Human
Expression cDNA Clone or AA Sequence:	Gly21 ^{Asn132}
Tag:	C-mFc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Note:	Recombinant Human Interleukin-13 is produced by our Mammalian expression system and the target gene encoding Gly21 ^{Asn132} is expressed with a mFc tag at the C-terminus.
Stability:	12 months from date of despatch
Summary:	Interleukin-13 is also known as IL-13. It is a protein that in humans is encoded by the IL13 gene. Interleukin-13 is an immunoregulatory cytokine produced primarily by activated Th2 cells.It is involved in several stages of B-cell maturation and differentiation. It up-regulates CD23 and MHC class II expression, and promotes IgE isotype switching of B cells. This cytokine down-regulates macrophage activity, thereby inhibits the production of pro-inflammatory cytokines and chemokines. This cytokine is found to be critical to the pathogenesis of allergen-induced asthma but operates through mechanisms independent of IgE and eosinophils.



[View online »](#)