

Product datasheet for TP726688

SLAMF7 Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant Human SLAMF7 (C-mFc)

Species: Human

Expression cDNA Clone

or AA Sequence:

Ser23-Met226

Tag: C-mFc

Buffer: Lyophilized from a 0.2 um filtered solution of PBS, pH7.4.

Note: Recombinant Human SLAM Family Member 7 is produced by our Mammalian expression

system and the target gene encoding Ser23-Met226 is expressed with a mFc tag at the C-

terminus.

Stability: 12 months from date of despatch

 Locus ID:
 57823

 UniProt ID:
 Q9NQ25

Summary: SLAMF7 is a single-pass type I membrane protein and contains 1 Ig-like C2-type

(immunoglobulin-like) domain. SLAMF7 is expressed in NK cells, activated B-cells, NK-cell line but not in promyelocytic, B-cell lines, or T-cell lines. Although the cytoplasmic domain of CS1 contains immunoreceptor tyrosine-based switch motifs (ITSM), which enables to recruite signaling lymphocyte activation molecule (SLAM)-associated protein (SAP/SH2D1A), it activates NK cells in the absence of a functional SAP. SLAMF7 positively regulated natural killer cell functions by a mechanism dependent on the adaptor EAT-2 but not the related adaptor SAP. However, in the absence of EAT-2, CRACC potently inhibited natural killer cell function. It was also inhibitory in T cells, which are typically devoid of EAT-2. Thus, SLAMF7 can exert activating or inhibitory influences on cells of the immune system depending on

cellular context and the availability of effector proteins.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com