

## Product datasheet for **TP726727**

### NOTCH2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Notch2 (C-6His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Leu26-Gln530
Tag:	C-6His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Note:	Recombinant Human Neurogenic Locus Notch Homolog Protein 2 is produced by our Mammalian expression system and the target gene encoding Leu26-Gln530 is expressed with a 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4853
UniProt ID:	<a href="#">Q04721</a>



[View online »](#)