

Product datasheet for **TP727010**

Cynomolgus Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Cynomolgus CD40 molecule/CD40(C-6His)
Species:	Cynomolgus
Expression cDNA Clone or AA Sequence:	Glu21-Arg193
Tag:	C-His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Note:	Recombinant Human CD40 molecule is produced by our expression system and the target gene encoding Glu21-Arg193 is expressed
Stability:	12 months from date of despatch
Summary:	CD40 is a Type I Transmembrane Glycoprotein that belongs to the TNF Receptor Superfamily. CD40 is expressed in B cells, follicular dendritic cells, dendritic cells, activated monocytes, macrophages, endothelial cells, vascular smooth muscle cells, and several tumor cell lines. The extracellular domain of CD40 is characterized by Cysteine rich repeat regions. Interaction of CD40 with its ligand (CD40L) leads to aggregation of CD40 molecules, which in turn interact with cytoplasmic components to initiate signaling pathways. Several different TRAF proteins (adaptor proteins) have been identified to serves as mediators of the signal transduction.



[View online »](#)