

Product datasheet for **TP727013**

Cynomolgus Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Cynomolgus Programmed cell death 1/PDCD1(C-Fc)
Species:	Cynomolgus
Expression cDNA Clone or AA Sequence:	Leu25-Gln167
Tag:	C-Fc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Note:	Recombinant Cynomolgus Programmed Cell Death Protein 1 is produced by our expression system and the target gene encoding Leu25-Gln167 is expressed
Stability:	12 months from date of despatch
Summary:	Programmed cell death protein 1(PDCD1) is a single-pass type I membrane protein and contains 1 Ig-like V-type domain. PD-1 is a member of the extended CD28/CTLA-4 family of T cell regulators. PDCD1 inhibits the T-cell proliferation and production of related cytokines including IL-1, IL-4, IL-10 and IFN- γ by suppressing the activation and transduction of PI3K/AKT pathway. In addition, coligation of PDCD1 inhibits BCR-mediating signal by dephosphorylating key signal transducer. PDCD1 has been suggested to be involved in lymphocyte clonal selection and peripheral tolerance, and thus contributes to the prevention of autoimmune diseases. As a cell surface molecule, PDCD1 regulates the adaptive immune response. Engagement of PD-1 by its ligands PD-L1 or PD-L2 transduces a signal that inhibits T-cell proliferation, cytokine production, and cytolytic function.



[View online »](#)