

Product datasheet for **TP727030**

IL13ra1 Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse IL-13RA1/NR4/CD213a1(C-AVI-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Ala25-Thr340
Tag:	C-AVI&His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4.
Note:	Biotinylated Recombinant Mouse Interleukin-13 receptor subunit alpha-1 is produced by our Mammalian expression system and the target gene encoding Ala25-Thr340 is expressed with a AVI, 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	16164
UniProt ID:	O09030
Summary:	Interleukin-13 receptor subunit alpha-1(IL-13RA1), also known as NR4, belongs to the hemopoietin receptor family. Interleukin-13 receptor is a complex of IL4R, IL13RA1, and possibly other components. It interacts with TRAF3IP1 and binds with low affinity to interleukin-13(IL3). IL-13RA1, together with IL4RA, can form a functional receptor for IL13. It also serves as an alternate accessory protein to the common cytokine receptor gamma chain for interleukin-4 (IL4) signaling, but cannot replace the function of IL2RG in allowing enhanced interleukin-2 (IL2) binding activity. The WSXWS motif in domains appears to be necessary for proper protein folding and thereby efficient intracellular transport and cell-surface receptor binding.



[View online »](#)