

## Product datasheet for TP727087

## **Mmp9 Mouse Recombinant Protein**

**Product data:** 

**Product Type: Recombinant Proteins** 

**Description:** Recombinant Mouse Matrix Metalloproteinase-9/MMP-9

Species: Mouse

**Expression cDNA Clone** 

or AA Sequence:

Ala20-Pro730

Tag: C-His

**Buffer:** Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.

Note: Recombinant Mouse Matrix Metalloproteinase-9 is produced by our Mammalian expression

system and the target gene encoding Ala20-Pro730 is expressed with a 10His tag at the C-

terminus.

Stability: 12 months from date of despatch

Locus ID: 17395 **UniProt ID:** P41245

Summary: Matrix metalloproteinases are a family of zinc and calcium dependent endopeptidases with

> the combined ability to degrade all the components of the extracellular matrix. MMP-9 (gelatinase B) can degrade a broad range of substrates including gelatin, collagen types IV and V, elastin and proteoglycan core protein. It is believed to act synergistically with interstitial collagenase (MMP1) in the degradation of fibrillar collagens as it degrades their denatured gelatin forms. MMP-9 is produced by keratinocytes, monocytes, macrophages and PMN leukocytes. MMP-9 is present in most cases of inflammatory responses. Structurally, MMP-9 may be divided into five distinct domains: a prodomain which is cleaved upon activation, a gelatinbinding domain consisting of three contiguous fibronectin type II units, a catalytic domain containing the zinc binding site, a prolinerich linker region, and a carboxyl

terminal hemopexinlike domain.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com