

Product datasheet for **TP727140**

Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse Host Cell Factor 2/HCFC2 (N-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Met1-Glu497
Tag:	N-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 2M Urea, pH 7.4.
Note:	Recombinant Mouse Host cell factor 2 is produced by our E.coli expression system and the target gene encoding Met1-Glu486 is expressed with a 6His tag at the N-terminus.
Stability:	12 months from date of despatch
Summary:	Host cell factor 2(HCFC2) is a cytoplasmic protein. It contains 2 fibronectin type-III domains.HCFC2 binds KMT2A/MLL1, as component of the MLL1/MLL complex.Hcfc2 negative regulation of transcription from RNA polymerase II promoter.



[View online »](#)