

Product datasheet for **TP727146**

ACYP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Acylphosphate Phosphohydrolase 1/ACYP1 (N-6His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Met1-Lys99
Tag:	N-His
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, 1mM DTT, pH 8.0.
Note:	Recombinant Human Acylphosphate phosphohydrolase 1 is produced by our E.coli expression system and the target gene encoding Met1-Lys99 is expressed with a 6His tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	97
UniProt ID:	P07311
Summary:	ACYP1, also known as Acylphosphatase-1, Acylphosphatase, erythrocyte isozyme, Acylphosphatase, organ-common type isozyme, Acylphosphate phosphohydrolase 1 and ACYPE, is a small cytosolic enzyme which catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates.ACYP1 is a protein which belongs to the acylphosphatase family and contains 1 fibrinogen C-terminal domain. Two isoenzymes have been isolated, called muscle acylphosphatase and erythrocyte acylphosphatase, on the basis of their tissue localization. This gene encodes the erythrocyte acylphosphatase isoenzyme. Alternatively spliced transcript variants that encode different proteins were identified through data analysis. Recombinant human ACYP1 protein was expressed in E. coli fused with HIS-tag at N-terminus.



[View online »](#)