

Product datasheet for **TP727288**

IL2 Receptor beta (IL2RB) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human IL-2 Receptor Subunit β /IL-2RB/CD122 (C-6His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Ala27-Asp239
Tag:	C-His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Note:	Recombinant Human Interleukin-2 Receptor subunit beta is produced by our Mammalian expression system and the target gene encoding Ala27-Asp239 is expressed with a 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3560
UniProt ID:	P14784
Summary:	Human IL-2RB, also known as interleukin-2 receptor subunit beta, is the receptor for interleukin-2. IL2 receptor complex is involved in receptor mediated endocytosis and transduces the mitogenic signals of IL2. IL2 receptor complex has three forms with respect to ability to bind IL2. IL-2RB is belonged to a type I membrane protein, and has a 26 residue signal peptide, a 214 residue extracellular region, a 25 residue transmembrane region and a 286 residue cytoplasmic domain. IL-2RB is the subunit critical for receptor-mediated signaling via physically or functionally coupling to other signaling molecules, such as the Jak-STAT and Src-family protein tyrosine kinase although it lacks apparent catalytic motifs.



[View online »](#)