

Product datasheet for **TP727291**

Car14 Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse Carbonic Anhydrase 14/CA14 (C-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Ala16-Met290
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Note:	Recombinant Mouse Carbonic anhydrase 14 is produced by our Mammalian expression system and the target gene encoding Ala16-Met290 is expressed with a 6His tag at the C-terminus.
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23831
UniProt ID:	Q9WVT6
Synonyms:	Carbonic Anhydrase 14; Carbonate Dehydratase XIV; Carbonic Anhydrase XIV; CA-XIV; CA14
Summary:	Mouse Ca14,also known as Carbonic anhydrase 14,is a member of large family of zinc metalloenzymes .It could catalyze reversible hydration of carbon dioxide. The reaction is fundamental to many processes such as respiration, renal tubular acidification and bone resorption. Fifteen CA isoforms have been reported so far. They have different patterns of tissue-specific expression and physiologic roles. Some CAs may serve as markers for tumors and hypoxia. CA XIV is a polypeptide consisting of an extracellular N-terminal catalytic domain, a membrane-spanning segment and a short intracellular C- terminal segment with several potential phosphorylation sites. A subset of CAs lack CA activity due to point mutations but retain esterase function. CA14 is widely expressed in the central nervous system



[View online »](#)