

Product datasheet for **TP727339**

Ces1d Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse Carboxylesterase 3/CES3 (C-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Tyr19-Glu561
Tag:	C-His
Buffer:	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Note:	Recombinant Mouse Carboxylesterase 3 is produced by our Mammalian expression system and the target gene encoding Tyr19-Glu561 is expressed with a 6His tag at the C-terminus.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	12 months from date of despatch
Locus ID:	104158
UniProt ID:	Q8VCT4
Summary:	Mouse Carboxylesterases 3 (CES3) is a member of five families of mammalian carboxylesterases that plays a role in catalyzing hydrolytic and transesterification reactions with xenobiotics, anticancer pro-drugs and narcotics, and detoxifying organophosphates and insecticides. Mammalian carboxylesterases are enzymes with broad substrate specificities ranging from small molecule esters to longchain fatty acid esters. It is shown that CESs has key roles in the metabolism of a wide variety of clinical drugs, illicit narcotics and chemical nerve agents. CES3 is broadly expressed in liver, colon and brain.



[View online »](#)