

Product datasheet for **TP727350**

Tnfrsf4 Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse OX40/TNFRSF4/CD134 (C-Fc)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Val20-Pro211
Tag:	C-Fc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Note:	Recombinant Mouse OX40L Receptor is produced by our Mammalian expression system and the target gene encoding Val20-Pro211 is expressed with a Fc tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	22163
UniProt ID:	P47741
Summary:	<p>OX40, also termed CD134 and TNFRSF4, is a T cell co-stimulatory molecule of the TNF receptor superfamily which plays a key role in the survival and homeostasis of effector and memory T cells. OX40 is expressed on CD4+ and CD8+ T cells upon engagement of the TCR by antigen presenting cells along with co-stimulation by CD40-CD40 Ligand and CD28-B7. The interaction between OX40 and OX40 ligand (OX40L) will occur when activated T cells bind to professional antigen-presenting cells (APCs). The T-cell functions, including cytokine production, expansion, and survival, are then enhanced by the OX40 costimulatory signals. OX40 signals are critical for controlling the function and differentiation of Foxp3+ regulatory T cells. OX40-OX40L interaction regulates T-cell tolerance, peripheral T-cell homeostasis, and T-cell-mediated inflammatory diseases.</p>



[View online »](#)