

## Product datasheet for **TP727403**

### Cd40 Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse CD40/TNFRSF5/CD40L Receptor (C-Fc)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Leu20-Arg193
Tag:	C-Fc
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Note:	Recombinant Mouse CD40L Receptor is produced by our Mammalian expression system and the target gene encoding Leu20-Arg193 is expressed with a Fc tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	21939
UniProt ID:	<a href="#">P27512</a>
Summary:	CD40 is a Type I Transmembrane Glycoprotein that belongs to the TNF Receptor Superfamily. CD40 is expressed in B cells, follicular dendritic cells, dendritic cells, activated monocytes, macrophages, endothelial cells, vascular smooth muscle cells, and several tumor cell lines. The extracellular domain of CD40 is characterized by Cysteine rich repeat regions. Interaction of CD40 with its ligand (CD40L) leads to aggregation of CD40 molecules, which in turn interact with cytoplasmic components to initiate signaling pathways. Several different TRAF proteins (adaptor proteins) have been identified to serves as mediators of the signal transduction. CD40 plays an essential role in mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development, and germinal center formation.



[View online »](#)