

## Product datasheet for **TP727452**

### Prlr Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse Prolactin Receptor/PRLR (C-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Gln20-Asp229
Tag:	C-His
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Note:	Recombinant Mouse Prolactin receptor is produced by our Mammalian expression system and the target gene encoding Gln20-Asp229 is expressed with a 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	19116
UniProt ID:	<a href="#">Q08501</a>
Summary:	The prolactin receptor (PRLR) is a member of the class I cytokine/lactogen receptor family which mediates the diverse cellular actions of prolactin in several tissues. PRLRs are expressed in normal and neoplastic human breast tissue, and in most breast cancer cells. PRLR contains an extracellular region that binds prolactin, a transmembrane region, and a cytoplasmatic region required for the activation of the Jak2-Stat5 signal transduction pathway by Prl which is essential for transcriptional activation of all known prolactin regulated genes. PRLRs have also been observed in ovarian follicular cells of mice, pigs, sheep, deer, and humans, as well as in luteal tissue in cow and horse ovaries. Furthermore, PRLR knockout mice exhibit failure of embryonic implantation, reduced number of mature oocytes, and low fertilization rates. Knockout females also display a reduced number of primary follicles.



[View online »](#)