

Product datasheet for **TP727708**

Slamf8 Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse SLAM Family Member 8/SLAMF8 (C-6His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	Val21-Asp231
Tag:	C-His
Buffer:	Lyophilized from a 0.2 um filtered solution of 20mMPB,150mMNaCl,pH7.4.
Note:	Recombinant Mouse SLAM Family Member 8 is produced by our Mammalian expression system and the target gene encoding Val21-Asp231 is expressed with a 6His tag at the C-terminus.
Storage:	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Stability:	12 months from date of despatch
Locus ID:	74748
UniProt ID:	Q9D3G2
Synonyms:	SLAM family member 8; B-lymphocyte activator macrophage expressed; CD353; Slamf8; Blame



[View online »](#)

Summary:

Mouse SLAM family member 8/SLAMF8 is a single-pass type I membrane protein. It contains one Ig-likeC2-type domain. SLAMF8 is expressed in lymphnode, spleen, thymus and bone marrow. The signaling lymphocyte activation molecule (SLAM) family includes homophilic and heterophilic receptors that modulate both adaptive and innate immune responses. These receptors share a common ectodomain organization: a membrane-proximal immunoglobulin constant domain and a membrane-distal immunoglobulin variable domain that is responsible for ligand recognition. SLAM family of receptors is expressed by a wide range of immune cells. Through their cytoplasmic domain, SLAM family receptors associate with SLAM-associated protein (SAP)-related molecules, a group of cytoplasmic adaptors composed almost exclusively of an SRC homology 2 domain. SLAM family receptors, in association with SAP family adaptors, have crucial roles during normal immune reactions in innate and adaptive immune cells. It may play a role in B-lineage commitment and/or modulation of signaling through the B-cell receptor.