

Product datasheet for **TP728070**

Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Biotinylated Human Nectin-2 (C-6His-Avi)
Species:	Human
Expression cDNA Clone or AA Sequence:	Gln32-Leu360
Tag:	C-6His-Avi
Buffer:	Lyophilized from a 0.2 um filtered solution of PBS, pH7.4.
Note:	Recombinant Human Poliovirus Receptor-Related Protein 2 is produced by our Mammalian expression system and the target gene encoding Gln32-Leu360 is expressed with a Avi, 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Summary:	CD112 is a type I transmembrane glycoprotein belonging to the Immunoglobulin superfamily. It comprises one Ig-like V-type domain and two Ig-like C2-type domains in the extracellular region. The V domain is believed to mediate nectin binding to its ligands. Nectin2 is known to bind the pseudorabies virus, and herpes simplex virus2 (HSV2), involving in cell to cell spreading of these viruses. It does not bind poliovirus. As a homophilic adhesion molecule, CD112 is found concentrated in adherens junctions, and exists on neurons, endothelial cells, epithelial cells and fibroblasts. CD112 has been identified as the ligand for DNAM-1 (CD226), and the interaction of CD226/CD112 mediates cytotoxicity and cytokine secretion by T and NK cells. The costimulatory responses may be a critical component in allergic reactions and may therefore become targets for anti-allergic therapy.



[View online »](#)