

Product datasheet for **TP728096**

SARS-CoV-2 Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant S1 Protein(Mammalian, C-6His)
Species:	SARS-CoV-2
Expression cDNA Clone or AA Sequence:	Val16-Arg685
Tag:	C-6His
Buffer:	Supplied as a 0.2 µm filtered solution of PBS, PH7.4.
Note:	Recombinant 2019-nCoV S1 Protein is produced by our Mammalian expression system and the target gene encoding Val16-Arg685 is expressed with a 6His tag at the C-terminus.
Stability:	12 months from date of despatch
Summary:	The spike (S) glycoprotein of coronaviruses is known to be essential in the binding of the virus to the host cell at the advent of the infection process. Most notable is severe acute respiratory syndrome (SARS). The severe acute respiratory syndrome-coronavirus (SARS-CoV) spike (S) glycoprotein alone can mediate the membrane fusion required for virus entry and cell fusion. It is also a major immunogen and a target for entry inhibitors. It's been reported that 2019-nCoV can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.



[View online »](#)