

## Product datasheet for **TP728143L**

### Recombinant Activin A, Human/Mouse/Rat

#### Product data:


<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Activin A, Human/Mouse/Rat
<b>Species:</b>	Human/Mouse/Rat
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	GLECDGKVNICCKKQFFVSKDIGWNDWIIAPSGYHANYCEGECPSHIAGTSGSSLSFHSTVINHYRMRGH SPFANLKSCCVPTKLRPMSMLYYDDGQNIKKDIQNMIVEECGCS
<b>Tag:</b>	Tag Free
<b>Predicted MW:</b>	The protein has a calculated MW of 13.1 kDa. The protein migrates as 15 kDa under reducing condition (SDS-PAGE analysis).
<b>Purity:</b>	>95% as determined by SDS-PAGE.
<b>Buffer:</b>	The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4.
<b>Bioactivity:</b>	Measure by its ability to inhibit the proliferation of mouse MPC-11 cells. The ED <sub>50</sub> for this effect is <10 ng/mL. The specific activity of recombinant human Activin A is approximately >1.0 x 10 <sup>3</sup> IU/mg.
<b>Endotoxin:</b>	<0.1 EU per 1 µg of the protein by the LAL method.
<b>Reconstitution Method:</b>	Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H <sub>2</sub> O to a concentration not less than 100 µg/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.
<b>Applications:</b>	Cell culture
<b>Storage:</b>	Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.
<b>UniProt ID:</b>	<a href="#">P08476</a>
<b>Synonyms:</b>	Inhibin beta A chain, Activin beta-A chain, Erythroid differentiation protein



[View online »](#)

**Summary:**

Activin is the member of the TGF beta superfamily of cytokines. The current understanding of Activin A is classified as a hormone, a Growth Factors, and a cytokine. Depending on the different cell type, that overexpression of Activin A can either inhibit or promote cells division. There are important implications for tumor biology. Activin A also increase joint inflammation in rheumatoid arthritis (RA) and induces pathogenic mechanisms in the respiratory system.

**Product images:**

SDS- PAGE analysis of recombinant human/mouse/rat Activin A