

Product datasheet for **TP728194L**

Recombinant FGF-14 (Fibroblast growth factor-14), Human

Product data:


| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant FGF-14 (Fibroblast growth factor-14), Human |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | AAAIASGLIRQKRQAREQHWDRPSASRRRSSPSKNGRLCNGNLVDIFS KVRIFGLKKRRLRRQDPQLKGIV TRLYCRQGYYLQMHPDGALDGTKDDSTNSTLFNLIPVGLRVVAIQGVKTGLYIAMNNGEGYLYPSELFTPEC KFKESVFENYYVIYSSMLYRQQESGRAWFLGLNKEGQAMKGNRVKTKPAAHFLPKPLEVAMYREPSLHD VGETVPKPGVTPSKSTSASAIMNNGGKPVNKSKT with polyhistidine tag at the N-terminus. |
| Tag: | His Tag (N-term) |
| Predicted MW: | The protein has a calculated MW of 28.28 kDa. The protein migrates as 33 kDa under reducing condition (SDS-PAGE analysis). |
| Purity: | >95% as determined by SDS-PAGE. |
| Buffer: | The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4. |
| Bioactivity: | Measure by its ability to induce 3T3 cells proliferation. The ED ₅₀ for this effect is <21 ng/mL. |
| Endotoxin: | <0.1 EU per 1 µg of the protein by the LAL method. |
| Reconstitution Method: | Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein. |
| Applications: | Cell culture |
| Storage: | Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles. |
| UniProt ID: | Q92915 |
| Synonyms: | FHF-4, FHF4, SCA27 |



[View online »](#)

Summary:

Fibroblast Growth Factors-14 (FGF-14) is a 27.7 kDa member of the fibroblast Growth Factors with 247 amino acid residues. FGF-14 is mainly expressed from brain, cervix. FGF-14 involved in nervous system development and function. May regulate voltage-gated sodium channels transport and function.

Product images:

SDS- PAGE analysis of recombinant human FGF-14