

Product datasheet for **TP728212M**

Recombinant Galectin-1, Human


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Galectin-1, Human
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	ACGLVASNLNLKPGECRLVRGEVAPDAKSFVLNLGKDSNNLCLHFNPRFNAHGDANTIVCNSKDGGAWGTEQREAVFPFQPGSVAEVCITFDQANLTVKLPDGYEFKFPNRLNLEAINYMAADGDFKIKCVAFD with polyhistidine tag at the N-terminus.
Tag:	His Tag (N-term)
Predicted MW:	The protein has a calculated MW of 15.5 kDa. The protein migrates as 15 kDa under reducing condition (SDS-PAGE analysis).
Purity:	>98% as determined by SDS-PAGE.
Buffer:	The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4.
Bioactivity:	Measured by its ability to agglutinate human red blood cells. The ED ₅₀ for this effect is <2 µg/mL.
Endotoxin:	<0.1 EU per 1 µg of the protein by the LAL method.
Reconstitution Method:	Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.
Applications:	Cell culture
Storage:	Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.
UniProt ID:	<u>P09382</u>
Synonyms:	Beta-galactoside-binding lectin L-14-I, Galaptin, 14 kDa lectin, S-LAC lectin-1


[View online »](#)

Summary:

Galectin-1 (Gal-1) is a lectin family member. Galectin-1 contains one carbohydrate recognition domain (CRD), responsible for β -galactoside binding, and is biologically active as homodimers. Several binding partners for galectin-1 have been identified, including CD45, CD43, CD69, pre-BCR, and vascular endothelial Growth Factors R2. Galectin-1 is constitutively presented in the immune system and other tissues, contributing to cell growth, proliferation, and apoptosis. The galectin-1 plays a central role in regulating acute and chronic inflammatory responses.

Product images:


SDS- PAGE analysis of recombinant human Galectin-1