

Product datasheet for **TP728268**

Recombinant IL-37 (Interleukin-37), Human

Product data:


| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant IL-37 (Interleukin-37), Human |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | MKNLNPKKFSIHDQDHKVLVLDSGNLI AVDPKNIYRPEIFFALASSLSSASAEKGSPI LLGVSKGEFCLYCDK DKGQSHPSLQLKKEKLMKLA AQKESARRPFIFYRAQVGSWNMLESA AHPGWFICTSCNCNEPVGVTDF ENRKHIEFSFQPVCKAEMSPSEVSD with polyhistidine tag at the C-terminus. |
| Tag: | His Tag (C-term) |
| Predicted MW: | The protein has a calculated MW of 19.49 kDa. The protein migrates as 22 kDa under reducing condition (SDS-PAGE analysis). |
| Purity: | >95% as determined by SDS-PAGE. |
| Buffer: | The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 8.0. |
| Bioactivity: | Measure by its ability to induce IL-8 secretion in human PBMCs. The ED ₅₀ for this effect is <0.9 µg/mL. |
| Endotoxin: | <0.1 EU per 1 µg of the protein by the LAL method. |
| Reconstitution Method: | Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein. |
| Applications: | Cell culture |
| Storage: | Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles. |
| UniProt ID: | Q9NZH6 |
| Synonyms: | IL-1F7, IL-1ζ (zeta), IL-1H4, FIL1, FIL1Z |



[View online »](#)

Summary:

Interleukin-37 (IL-37) consists of 192 amino acids. IL-37 is a member of the interleukin-1 cytokine family. The primary function is binding to the interleukin-18 receptor (IL18R1 /IL-1Rrp) and then becoming a subunit to inhibit the activity of IL-18. This behavior often occurs between various immune responses. It also significantly affects the negative regulation of tumor necrosis factor production.

Product images:

SDS- PAGE analysis of recombinant human IL-37