

Product datasheet for **TP728290**

Recombinant TNF beta (Tumor necrosis factor beta), Human

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant TNF beta (Tumor necrosis factor beta), Human
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	LPGVGLTPSAAQTARQHHPKMHLAHLSTLKPAAHLIGDPSKQNSLLWRANTDRAFLQDGFSLSNNSLLVPTSGIYFVYSQVWFSGKAYSPKATSSPLYLAHEVQLFSSQYPFHVPLLSSQKMVYPGLQEPWLHSMYHGAAFAQLTQGDQLSTHTDGIPIHLVLSPTVFFGAFAL with polyhistidine tag at the N-terminus.
Tag:	His Tag (N-term)
Predicted MW:	The protein has a calculated MW of 18.65 kDa. The protein migrates as 21 kDa under reducing condition (SDS-PAGE analysis).
Purity:	>98% as determined by SDS-PAGE.
Buffer:	The protein was lyophilized from a 0.2 µm filtered solution containing 20 mM sodium citrate, 0.2 M NaCl, pH 3.5.
Bioactivity:	Measure by its ability to induce cytotoxicity in L929 cells in the presence of actinomycin D. The ED ₅₀ for this effect is <3 pg/mL. The specific activity of recombinant human TNF beta is > 3.3x 10 ⁸ IU/mg.
Endotoxin:	<0.1 EU per 1 µg of the protein by the LAL method.
Reconstitution Method:	Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.
Applications:	Cell culture
Storage:	Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.
UniProt ID:	<u>P01374</u>
Synonyms:	TNFSF1B, Lymphotoxin-alpha (LT-α),LTA



[View online »](#)

Summary:

Human TNF beta (Lymphotoxin-alpha) is a 22 kDa cytokine with 205 amino acid residues. TNF beta belongs to tumor necrosis factor (TNF) superfamily as TNF alpha. Binding to lymphotoxin- β receptor (Lt β R), canonical and non-canonical NF- κ B signal pathway will be activated. TNF beta is related to apoptosis and inflammatory signals

Product images:

SDS- PAGE analysis of recombinant human TNF beta