


## Product datasheet for TP750180

### Recombinant SARS coronavirus 2 Nucleocapsid Protein (N-protein)

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant SARS coronavirus 2 Nucleocapsid Protein (N-protein) with N-terminal HIS tag, expressed in E.Coli
Species:	SARS-CoV-2
Expression Host:	E. coli
Tag:	N-HIS
Predicted MW:	47 kDa
Concentration:	>50 ug/mL as determined by Bradford protein assay method.
Purity:	> 90% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25mM Tris, pH8.0, 250mM NaCl, 10% glycerol.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">QHD43423</a></u>

#### Product images:



Purified recombinant protein N protein was analyzed by SDS-PAGE gel and Coomassie Blue Staining.



[View online »](#)