


Product datasheet for TP750182

SARS-CoV-2 Spike-RBD Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant SARS coronavirus 2 Spike-RBD protein(Lys310-Leu560), with N-terminal His tag, expressed in E.coli cells, 50ug
Species:	SARS-CoV-2
Expression Host:	E. coli
Tag:	N-His
Predicted MW:	28.1 kDa
Concentration:	>50 ug/mL as determined by Bradford Protein Assay method
Purity:	> 90% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50mM Tris,pH8.0,500mM NaCl,8M Urea.
Storage:	Store at +4°C.
Stability:	Stable for at least 1 year from receipt of products under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>QHD43416</u>

Product images:



Purified recombinant protein S Protein (RBD) was analyzed by SDS-PAGE gel and Coomassie Blue Staining.



[View online »](#)