

Product datasheet for **TP760165**

HDGFL3 (NM_016073) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hepatoma-derived growth factor, related protein 3 (HDGFRP3), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length HDGFRP3
Tag:	N-His
Predicted MW:	22.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057157
Locus ID:	50810
UniProt ID:	Q9Y3E1 , A0A024R216
RefSeq Size:	1973
Cytogenetics:	15q25.2
RefSeq ORF:	609
Synonyms:	CGI-142; HDGF-2; HDGF2; HDGFRP3; HRP-3
Summary:	Enhances DNA synthesis and may play a role in cell proliferation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images: