

## Product datasheet for **TP760179**

### **C19orf48 (NM\_199249) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human chromosome 19 open reading frame 48 (C19orf48), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence encoding human full-length C19orf48
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	12.9 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_954857</a>
<b>Locus ID:</b>	84798
<b>UniProt ID:</b>	<a href="#">Q6RUI8</a> , <a href="#">A0A024R4G9</a>
<b>RefSeq Size:</b>	1828
<b>Cytogenetics:</b>	19q13.33
<b>RefSeq ORF:</b>	351



[View online »](#)

**Product images:**

