

## **Product datasheet for TP760196**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## PRPK (TP53RK) (NM 033550) Human Recombinant Protein

**Product data:** 

**Product Type: Recombinant Proteins** 

Description: Recombinant protein of human TP53 regulating kinase (TP53RK), full length, with N-terminal

HIS tag, expressed in E.Coli, 50ug

Human Species: **Expression Host:** E. coli

**Expression cDNA Clone** 

or AA Sequence:

A DNA sequence encoding human full-length TP53RK

N-His Tag:

Predicted MW: 28 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining **Buffer:** 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Store at -80°C. Storage:

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 291028

Locus ID: 112858 UniProt ID: Q96S44 3384 RefSeq Size: 20q13.12

Cytogenetics:

RefSeq ORF: 759

Synonyms: BUD32; C20orf64; dJ101A2; GAMOS4; Nori-2; Nori-2p; PRPK; TPRKB





**Summary:** 

Component of the EKC/KEOPS complex that is required for the formation of a threonylcarbamoyl group on adenosine at position 37 (t(6)A37) in tRNAs that read codons beginning with adenine (PubMed:22912744, PubMed:27903914). The complex is probably involved in the transfer of the threonylcarbamoyl moiety of threonylcarbamoyl-AMP (TC-AMP) to the N6 group of A37 (PubMed:22912744, PubMed:27903914). TP53RK has ATPase activity in the context of the EKC/KEOPS complex and likely plays a supporting role to the catalytic subunit OSGEP (By similarity). Atypical protein kinase that phosphorylates 'Ser-15' of p53/TP53 protein and may therefore participate in its activation (PubMed:11546806). [UniProtKB/Swiss-Prot Function]

**Protein Families:** 

Druggable Genome, Protein Kinase

## **Product images:**

