

Product datasheet for TP760201

Cytokeratin 18 (KRT18) (NM_199187) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human keratin 18 (KRT18), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Tag:	N-His
Predicted MW:	47.9 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_954657
Locus ID:	3875
RefSeq Size:	1439
Cytogenetics:	12q13.13
RefSeq ORF:	1290
Synonyms:	CK-18; CYK18; K18
Summary:	KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
Protein Pathways:	Pathogenic Escherichia coli infection



[View online »](#)

Product images:

