

Product datasheet for **TP760208**

HMGB4 (NM_001008728) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human high-mobility group box 4 (HMGB4), transcript variant 2, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length HMGB4
Tag:	N-His
Predicted MW:	13.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001008728
Locus ID:	127540
UniProt ID:	Q8WW32
RefSeq Size:	637
Cytogenetics:	1p35.1
RefSeq ORF:	336
Synonyms:	dj1007G16.5; FLJ40388; MGC88128
Protein Families:	Transcription Factors



[View online »](#)

Product images: