

Product datasheet for **TP760302**

ZNF354A (NM_005649) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 354A (ZNF354A), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length ZNF354A
Tag:	N-His
Predicted MW:	69.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005640
Locus ID:	6940
UniProt ID:	O60765 , V9HWI2
RefSeq Size:	2555
Cytogenetics:	5q35.3
RefSeq ORF:	1815
Synonyms:	EZNF; HEL104; HKL1; KID-1; KID1; TCF17
Protein Families:	Transcription Factors



[View online »](#)

Product images: