

Product datasheet for **TP760485**

MFAP1 (NM_005926) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human microfibrillar-associated protein 1 (MFAP1), with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length MFAP1
Tag:	N-His
Predicted MW:	51.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005917
Locus ID:	4236
UniProt ID:	P55081
RefSeq Size:	2063
Cytogenetics:	15q15.3
RefSeq ORF:	1317
Synonyms:	AMP
Summary:	Involved in pre-mRNA splicing as a component of the spliceosome.[UniProtKB/Swiss-Prot Function]
Protein Families:	Secreted Protein



[View online »](#)

Product images:

