

## Product datasheet for **TP760509**

### ZNF596 (NM\_001042415) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human zinc finger protein 596 (ZNF596), transcript variant 2, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length ZNF596
Tag:	N-His
Predicted MW:	58.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001035880</a>
Locus ID:	169270
UniProt ID:	<a href="#">Q8TC21</a>
RefSeq Size:	2730
Cytogenetics:	8p23.3
RefSeq ORF:	1512
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**Product images:**