

Product datasheet for **TP760517**

CCNB1IP1 (NM_021178) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Human cyclin B1 interacting protein 1, E3 ubiquitin protein ligase (CCNB1IP1), transcript variant 1, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | A DNA sequence encoding human full-length CCNB1IP1 |
| Tag: | N-His |
| Predicted MW: | 31.4 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1% sarkosyl, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_067001 |
| Locus ID: | 57820 |
| UniProt ID: | Q9NPC3 , A0A384MR52 |
| RefSeq Size: | 1499 |
| Cytogenetics: | 14q11.2 |
| RefSeq ORF: | 831 |
| Synonyms: | C14orf18; HEI10 |
| Summary: | HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.[supplied by OMIM, Apr 2004] |



[View online »](#)

Product images: