

## Product datasheet for **TP760540**

### PRELID2 (NM\_138492) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human PRELI domain containing 2 (PRELID2), transcript variant 3, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length PRELID2
Tag:	N-His
Predicted MW:	17 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_612501</a>
Locus ID:	153768
UniProt ID:	<a href="#">Q8N945</a>
RefSeq Size:	5047
Cytogenetics:	5q32
RefSeq ORF:	444



[View online »](#)

**Product images:**

