

Product datasheet for **TP760560**

DHRS13 (NM_144683) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human dehydrogenase/reductase (SDR family) member 13 (DHRS13), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	A DNA sequence encoding human full-length DHRS13
Tag:	N-His
Predicted MW:	40.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, pH 8.0, 8 M urea
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_653284
Locus ID:	147015
UniProt ID:	Q6UX07
RefSeq Size:	1939
Cytogenetics:	17q11.2
RefSeq ORF:	981
Synonyms:	SDR7C5
Summary:	Putative oxidoreductase.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome


[View online »](#)

Product images:

