

Product datasheet for **TP760596**

H3.3A (H3F3A) (NM_002107) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Human H3 histone, family 3A (H3F3A), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | A DNA sequence encoding human full-length H3F3A |
| Tag: | N-His |
| Predicted MW: | 15.1 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 50 mM Tris-HCl, pH 8.0, 8 M urea |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_002098 |
| Locus ID: | 3020 |
| UniProt ID: | P84243 , B2R4P9 |
| RefSeq Size: | 2263 |
| Cytogenetics: | 1q42.12 |
| RefSeq ORF: | 408 |
| Synonyms: | H3-3B; H3.3A; H3F3; H3F3A |



[View online »](#)

Summary:

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene contains introns and its mRNA is polyadenylated, unlike most histone genes. The protein encoded is a replication-independent member of the histone H3 family. [provided by RefSeq, Jul 2008]

Protein Pathways:

Systemic lupus erythematosus

Product images: